

PRMT9 Rabbit pAb

货号: B20078

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:1000 - 1:2000
理论分子量	81kDa/94kDa
实测分子量	105kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,HeLa,SGC-7901
细胞定位	cytoplasm
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-170 of human PRMT9 (NP_612373.2).
序列	MSNSRPRSRRDAGGGAGAAGRDELVSRLQSAEHCLGVQDFGTAYAHYLLVLSLAPELKHDKVETFYTLFRWAEELDA LSRIQDLLGCYEQALELFPDDEVICNSMGEHLFRMGFRDEAAGYFHKAVKLNPDFSDAKENFYRVANWLVERWHFIMLND TKRNTIYNAAI

靶点信息

研究背景	This gene encodes a type II methyltransferase. Post-translational modification of target proteins by PRMT s plays an important regulatory role in many biological processes, whereby PRMTs methylate arginine residues by transferring methyl groups from S-adenosyl-L-methionine to the guanidino nitrogen atoms of arginine. The protein encoded by this gene methylates spliceosome associated protein 145 to regulate alternative splicing and acts as a modulator of small nuclear ribonucleoprotein maturation. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.
------	--

基因ID	90826
基因名	PRMT9
Swiss	Q6P2P2
别名	PRMT9;PRMT10

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn