

DNAJC5 Rabbit pAb

货号: B20080

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:1000 - 1:2000 IHC: 1:50 - 1:200
理论分子量	18kDa/22kDa
实测分子量	30kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain,Mouse liver,Rat brain
细胞定位	Cell membrane,Lipid-anchor,Melanosome,Membrane
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-198 of human DNAJC5 (NP_079495.1).
序列	MADQRQRSLSTSGESLYHVLGLDKNATSDDIKSYRKLALKYHPDKNPDNPEAADKFKEINNAHAILTDATKRNIYDKYGS LGLYVAEQFGEENVNTYFVLSSWWAKALFVFCGLLTCCYCCCCCLCCCFNCCCGKCKPKAPEGEETEFYVSPEDLEAQLQS DEREATDTPIVIQPASATETTQLTADSHPSYHTDGFN

靶点信息

研究背景	This gene is a member of the J protein family. J proteins function in many cellular processes by regulating the ATPase activity of 70 kDa heat shock proteins. The encoded protein plays a role in membrane trafficking and protein folding, and has been shown to have anti-neurodegenerative properties. The encoded protein is known to play a role in cystic fibrosis and Huntington's disease. A pseudogene of this gene is located on the short arm of chromosome 8.
------	--

基因ID	80331
基因名	DNAJC5
Swiss	Q9H3Z4
别名	DNAJC5;CLN4;CLN4B;CSP;DNAJC5A;NCL;mir-941-2;mir-941-3;mir-941-4;mir-941-5

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn