

# ATF7IP Rabbit pAb

货号: B20092

## 产品信息

|       |                                                                                                        |
|-------|--------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse                                                                                            |
| 宿主    | Rabbit                                                                                                 |
| 克隆性   | Polyclonal                                                                                             |
| 预测反应  |                                                                                                        |
| 应用    | WB                                                                                                     |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:1000                                                                              |
| 理论分子量 | 118kDa/136kDa/137kDa                                                                                   |
| 实测分子量 | 250KDa                                                                                                 |
| 形式    | Liquid                                                                                                 |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                           |
| 阳性对照  | MCF7,A-549,Mouse kidney                                                                                |
| 细胞定位  | Nucleus                                                                                                |
| 纯化    | Affinity purification                                                                                  |

## 抗原信息


|      |                                                                                                                                                                                                                 |
|------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 525-725 of human ATF7 IP (NP_060649.3).                                                                                           |
| 序列   | NEKDNKPEEEEQVIHEDDERPSEKNEFSRRKRKSKSEMDNVQSKRRRYMEEYEAEFQVKITAKGDINQKLQKVIQWLLE<br>EKLALQCAVFDKTLAELKTRVEKIECNKRHKTVLTELQAKIARLTKRFEAAKEDLKKRHEHPPNPPVSPGKTVNDVNSNNN<br>MSYRNAGTVRQMLESKRNVSESAPPSFQTPVNTVSSTNL |

## 靶点信息

|      |                                                                                                                                                                                                                                                               |
|------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | ATF7IP is a multifunctional nuclear protein that associates with heterochromatin. It can act as a transcriptional coactivator or corepressor depending upon its binding partners (summary by Liu et al., 2009 [PubMed 19106100]).[supplied by OMIM, Nov 2010] |
| 基因ID | 55729                                                                                                                                                                                                                                                         |

|       |                                  |
|-------|----------------------------------|
| 基因名   | ATF7IP                           |
| Swiss | Q6VMQ6                           |
| 别名    | ATF7IP;AM;ATF-IP;MCAF;MCAF1;p621 |

产品验证



Western blot analysis of ATF7IP expressed in MCF7,A-549,Mouse kidney using ATF7IP Rabbit pAb at 1:1 000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. B locking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)