

CHFR Rabbit pAb

货号: B20093

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	63kDa/69kDa/72kDa/73kDa
实测分子量	70kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HT-1080,22Rv1,HepG2,Mouse kidney,Mouse testis
细胞定位	Nucleus,PML body
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-290 of human CHFR (NP_001154818.1).
序列	MERPEEGKQSPPPQPWGRLLRLGAEEGEPHVLLRKREWTIGRRRGCDLSFPSNKLVS GDHCRIVVDEKSGQVTLEDTSTS GTVINKLKVVKKQTCPLQTGDVIYLVYRKNEPEHNVAYLYESLSEKQGMTQESFDTSGAGAGRGADPRVPPSSPATQVCF EEPQPSTSTSDLFPTASASSTEPSPAGRERSSSCGSGGGGISPKGSGPSVASDEVSSFASALPDRKTASFSSLEPQDQEDL EPVKKKMRGDGDLNLNGQLLVAQPRRNAQTVHEDVRAAAGKPKMEETL

靶点信息

研究背景	This gene encodes an E3 ubiquitin-protein ligase required for the maintenance of the antephase checkpoint that regulates cell cycle entry into mitosis and, therefore, may play a key role in cell cycle progression and tumorigenesis. The encoded protein has an N-terminal forkhead-associated domain, a central RING-finger domain, and a cysteine-rich C-terminal region. Alternatively spliced transcript variants that encode different protein isoforms have been described.
------	--

基因ID	55743
基因名	CHFR
Swiss	Q96EP1
别名	CHFR;RNF116;RNF196

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn