

GNA14 Rabbit pAb

货号: B20094

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	41kDa
实测分子量	41kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	BxPC-3,293T,Mouse stomach,Mouse lung,Rat kidney
细胞定位	extracellular exosome,plasma membrane
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 40-350 of human GNA14 (NP_004288.1).
序列	LLGTGESGKSTFIKQMRIIHGSGYSDEDRKGFTKL VYQNIFTAMQAMIRAMDTLRIQYVCEQNKENAQIIREVEVDKVSMLSREQVEAIKQLWQDPG IQECYDRRREYQLSDSAKYLT DIDRIATPSFVPTQQDVL RVRVPTTGII EYPFDLENIIFRMVDVGGQRSERRKWIHCFESVTSIIFLVALSEYDQVLAECDNENRMEESKALFKTIITYPWFLNSSVILFLNKKDLL EEKIMYSHLISYFPEYTGPKQDVRAARDFILKLYQDQNP DKEKVIYSHFTCATDTDNIRFVFAAVKDTILQLNLR

靶点信息

研究背景	This gene encodes a member of the guanine nucleotide-binding, or G protein family. G proteins are heterotrimers consisting of alpha, beta and gamma subunits. The encoded protein is a member of the alpha family of G proteins, more specifically the alpha q subfamily of G proteins. The encoded protein may play a role in pertussis-toxin resistant activation of phospholipase C-beta and its downstream effectors.
------	---

基因ID	9630
基因名	GNA14
Swiss	O95837
别名	GNA14

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn