

RAB34 Rabbit pAb

货号: B20116

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	21kDa/28kDa/29kDa
实测分子量	29kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse lung,Mouse brain,Mouse kidney,Mouse testis
细胞定位	Nucleus,nucleolus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-259 of human RAB34 (NP_114140.4).
序列	MNILAPVRRDRVLAELPQCLRKEAALHGHKDFHPRVTCACQEHRTGTVGFKISKVIVVGDL SVGKTCLINRFCKDTFDKNY KATIGVDFEMERFEVLGIPFSLQLWDTAGQERFKCIASTYYRGAQAIIVFNLNDVASLEHTKQWLADALKENDPSSVLLFLV GSKKDLSTPAQYALMEKDALQVAQEMKAEYWAVSSLTGENVREFFRVAALTFEANVLAELEKSGARRIGDVVRINSDDS NLYLTASKKKPTCCP

靶点信息

研究背景	This gene encodes a protein belonging to the RAB family of proteins, which are small GTPases involved in protein transport. This family member is a Golgi-bound member of the secretory pathway that is involved in the repositioning of lysosomes and the activation of macropinocytosis. Alternative splicing of this gene results in multiple transcript variants. An alternatively spliced transcript variant produces the nine-amino acid residue-repeats (NARR) protein, which is a functionally distinct nucleolar protein resulting from a different reading frame.
------	---

基因ID	83871
基因名	RAB34
Swiss	P0DI83,Q9BZG1
别名	RAB34;NARR;RAB39;RAH

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn