

L3MBTL2 Rabbit pAb

货号: B20119

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:1000 - 1:2000
理论分子量	69kDa/79kDa
实测分子量	95kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,HL-60,HT-29
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 486-705 of human L3M BTL2 (NP_113676.2).
序列	PATFCQKNDIELTPPKGYESAQTFNWENYLEKTKSKAAPSRLFNMDCPNHGFKVGMKLEAVDLMPEPRLICVATVKRVVHRLLSIHFDGWDSEYDQWVDCESPDIYPVGWCELTGYQLQPPVAAEPATPLKAKEATKKKKKQFGKKRKRIPTKTRPLRQGSKKPLLEDDPQGARKISSEPVPGEIHAVRVKEEHLDVASPDKASSPELPVSVENIKQETDD

靶点信息

研究背景	Putative Polycomb group (PcG) protein. PcG proteins maintain the transcriptionally repressive state of genes, probably via a modification of chromatin, rendering it heritably changed in its expressibility. Its association with a chromatin-remodeling complex suggests that it may contribute to prevent expression of genes that trigger the cell into mitosis. Binds to monomethylated and dimethylated 'Lys-20' on histone H4. Binds histone H3 peptides that are monomethylated or dimethylated on 'Lys-4', 'Lys-9' or 'Lys-27'.
------	--

基因ID	83746
基因名	L3MBTL2
Swiss	Q969R5
别名	L3MBTL2;H-I(3)mbt-I;L3MBT

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn