

ENC1 Rabbit pAb

货号: B20151

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:1000 - 1:2000
理论分子量	57kDa/66kDa
实测分子量	66kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,U-87MG,SW480,MCF7,HT-1080,Mouse brain,Rat brain
细胞定位	Cytoplasm,Nucleus matrix,cytoskeleton
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human ENC1 (NP_003624.1).
序列	MSVSVHENRKSRASSGSINIYLFHKSSYADSVLTHLNLLRQQRLFTDVLLHAGNRTFPCHRAVLAACSRYEAMFSGGLKE SQDSEVNFDSNIHPEVLELLLDYAYSSRVIINEENAESLLEAGDMLEFQDIRDACAEFLEKNLHPTNCLGMLLLSDAHQCTKL YELSWRMCLSNFQTIRKNEDFLQLPQDMVVQLLSSEELETEDERLVYESAINWISYDLKKRYCYLPELLQTVRLALLPAIYLM ENVAMEELITKQRKSKEIVEEAI

靶点信息

研究背景	This gene encodes a member of the kelch-related family of actin-binding proteins. The encoded protein pl ays a role in the oxidative stress response as a regulator of the transcription factor Nrf2, and expression o f this gene may play a role in malignant transformation. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.
------	---

基因ID	8507
基因名	ENC1
Swiss	O14682
别名	ENC1;CCL28;ENC-1;KLHL35;KLHL37;NRPB;PIG10;TP53I10

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn