

SMARCA1 Rabbit pAb

货号: B20154

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:1000 - 1:2000 IHC: 1:50 - 1:200
理论分子量	121kDa/122kDa
实测分子量	140kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,293T,A-549,HT-1080,LO2
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human SMARCA1 (NP_003060.2).
序列	MEQDTAAVAATVAAADATATIVVIEDEQPGPSTSQEEGAAAAATEATAATEKGEKKKEKNVSSFQLKLAAPKSEKEMDPEYEEKMKADRAKRFEFLKQTELFHFQPSAQKSPTSPLNMKLGRPRIKKDEKQSLISAGDYRHRRTQEEDDELLSESRTKSNVCIRFEVSPSYVKGGLRDYQIRGLNWLISLYEN

靶点信息

研究背景	This gene encodes a member of the SWI/SNF family of proteins. The encoded protein is an ATPase which is expressed in diverse tissues and contributes to the chromatin remodeling complex that is involved in transcription. The protein may also play a role in DNA damage, growth inhibition and apoptosis of cancer cells. Alternative splicing results in multiple transcript variants.
------	--

基因ID	6594
基因名	SMARCA1
Swiss	P28370
别名	SMARCA1;ISWI;NURF140;SNF2L;SNF2L1;SNF2LB;SNF2LT;SWI;SWI2;hSNF2L

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn