

ZSCAN4C Rabbit pAb

货号: B20168

产品信息

反应	Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	57kDa
实测分子量	57kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse kidney
细胞定位	Chromosome,Nucleus,telomere
纯化	Affinity purification

抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 400-506 of mouse ZSCAN4 (NP_001013787.1).
序列	CSRMFKHARSLSSHQRTLNNKSELLCVTCQKMFKRVSDDRTHEIIHMPEKPFKCSTCEKSFSHKTNLKSHEMIHTGEMPYVCSLCSRRFRQSSTYHRHLRNYHRSD

靶点信息

研究背景	Predicted to enable DNA-binding transcription factor activity, RNA polymerase II-specific and RNA polymerase II cis-regulatory region sequence-specific DNA binding activity. Acts upstream of or within negative regulation of mitotic recombination; regulation of transcription by RNA polymerase II; and telomere maintenance via telomere lengthening. Located in chromosome, telomeric region; cytoplasm; and nucleus. Is expressed in 1-cell stage embryo; 2-cell stage embryo; primary oocyte; and secondary oocyte. Orthologous to human ZSCAN4 (zinc finger and SCAN domain containing 4).
------	--

基因ID	245109
基因名	
Swiss	
别名	ZSCAN4;ZNF494;ZSCAN4C

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn