

IL17RB Rabbit pAb

货号: B20186

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:1000 - 1:4000
理论分子量	31kDa/55kDa
实测分子量	75kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,293T,SW620,A-549,B-cell,Mouse liver,Rat liver
细胞定位	Cell membrane,Secreted,Single-pass type I membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 18-292 of human IL17R B (NP_061195.2).
序列	REPTVQCGSETGPSPEWMLQHDLPGLDLRDLRVEPVTTSVATGDYSILMNVSWVLRADASIRLLKATKICVTGKSNFQSYS CVRRCNYTEAFQTQTRPSGGKWTFSYIGFPVELNTVYFIGAHNIPNANMNEDGPSMSVNFTSPGCLDHIMKYKKKCVKAGS LWDPNITACKKNEETVEVNFTTTPLGNRYMALIQHSTIIGFSQVFEPHQKKQTRASVVIPVTGDSEGATVQLTPYFPTCGSD CIRHKGTVVLCPCQTGVFPFLDNNKSKPGGWLP

靶点信息

研究背景	The protein encoded by this gene is a cytokine receptor. This receptor specifically binds to IL17B and IL17 E, but does not bind to IL17 and IL17C. This receptor has been shown to mediate the activation of NF-kap paB and the production of IL8 induced by IL17E. The expression of the rat counterpart of this gene was fo und to be significantly up-regulated during intestinal inflammation, which suggested the immunoregulatory activity of this receptor.
------	---

基因ID	55540
基因名	IL17RB
Swiss	Q9NRM6
别名	IL17RB;CRL4;EVI27;IL17BR;IL17RH1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn