

MESP1 Rabbit pAb

货号: B20187

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:200 - 1:1000
理论分子量	28kDa
实测分子量	20kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,HeLa,THP-1,Mouse heart,Rat heart
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 149-268 of human MES P1 (NP_061140.1).
序列	RRCRQRGDAGSPRGCLPCDDCPAQMQRTRQAEGQGQGRGLGLVSAVRAGASWGSPACPGARAAPEPRDPPALFAE AACPEGQAMEPSPPSPLLPGDVLALLETWMPLSPLEWLPEEPK

靶点信息

研究背景	Transcription factor. Plays a role in the epithelialization of somitic mesoderm and in the development of c ardiac mesoderm. Defines the rostrocaudal patterning of the somites by participating in distinct Notch pat hways (By similarity.
基因ID	55897

基因名	MESP1
Swiss	Q9BRJ9
别名	MESP1;bHLHc5

产品验证



Western blot analysis of MESP1 expressed in Mouse heart,Rat heart,A-549 using MESP1 Rabbit pAb at 1 :1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane . Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn