

CECR2 Rabbit pAb

货号: B20190

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:200 - 1:1000
理论分子量	58kDa/162kDa/164kDa
实测分子量	164kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse liver,Rat liver
细胞定位	nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 170-370 of human CECR2 (NP_113601.2).
序列	VQGKSNGELSLSRESEGQKNVSSIPGKTGKRRGRPPKRKKLQEEILLSEKQEENSLASEPQTRHGSQGPQGQGTWWLLCQ TEEEWRQVTESFRERTSLRERQLYKLLSEDFLPEICNMIAQKGKRPQRTKAELHPRWMSDHLSIKPVKQEETPVLTRIEKQK RKEEEEERQILLAVQKKEQEQLKEERKRELEEKVKAVEG

靶点信息

研究背景	This gene encodes a bromodomain-containing protein that is involved in chromatin remodeling, and may additionally play a role in DNA damage response. The encoded protein functions as part of an ATP-dependent complex that is involved in neurulation. This gene is a candidate gene for Cat Eye Syndrome. Alternative splicing results in multiple transcript variants.
基因ID	27443

基因名	CECR2
Swiss	Q9BXF3
别名	CECR2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn