

GABRR2 Rabbit pAb

货号: B20206

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	54kDa/56kDa
实测分子量	65kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HL-60,HepG2,Mouse heart
细胞定位	Cell junction,Cell membrane,Multi-pass membrane protein,postsynaptic cell membrane,synapse
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 350-450 of human GABRR2 (NP_002034.3).
序列	TVQERKERKLREKFPCMCGMLHSKTMMLDGSYSESEANSLAGYPRSHILTEEERQDKIVVHLGLSGEANAARKKGLLKGQ TGFRIFQNTHAIDKYSRLIFP

靶点信息

研究背景	Gamma-aminobutyric acid (GABA) is the major inhibitory neurotransmitter in the mammalian brain where it acts at GABA receptors, which are ligand-gated chloride channels. The protein encoded by this gene is a member of the rho subunit family and is a component of the GABA type A receptor complex. This gene exists on chromosome 6q next to the gene encoding the rho 1 subunit of the GABA type A receptor, in a region thought to be associated with susceptibility for psychiatric disorders and epilepsy. Polymorphisms in this gene may also be associated with alcohol dependence, and general cognitive ability.
------	--

基因ID	2570
基因名	GABRR2
Swiss	P28476
别名	GABRR2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn