

HGFAC Rabbit pAb

货号: B20208

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	70kDa
实测分子量	75kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	LO2,BT-474,SW620,Mouse liver,Mouse pancreas,Rat liver
细胞定位	Secreted
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 408-655 of human HGF AC (NP_001519.1).
序列	IIGGSSSLPGSHPWLAAYIGDSFCAGSLVHTCWVVSAAHCFSHSPPRDSVSVVLGQHFFNRTTDTQTFGIEKYIPYTLYS VFNPSDHDLVLIRLKKKGDRCATRSQFVQPICLPEPGSTFPAGHKCQIAGWGHLDENVSGYSSSLREALVPLVADHKCSS PEVYGADISPNNMLCAGYFDCKSDACQGDSSGGPLACEKNGVAYLYGIISWGDGCGRLHKPGVYTRVANYVDWINDRIRPP RRLVAPS

靶点信息

研究背景	This gene encodes a member of the peptidase S1 protein family. The encoded protein is first synthesized as an inactive single-chain precursor before being activated to a heterodimeric form by endoproteolytic processing. It acts as serine protease that converts hepatocyte growth factor to the active form. Alternatively splicing results in multiple transcript variants.
------	---

基因ID	3083
基因名	HGFAC
Swiss	Q04756
别名	HGFAC;HGFA

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn