

ITGA11 Rabbit pAb

货号: B20218

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	133kDa
实测分子量	160kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,SH-SY5Y
细胞定位	Membrane,Single-pass type I membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 940-1141 of human ITG A11 (NP_001004439.1).
序列	KEDNVAPLRFHCLKYEADVLFTRSSSLSHYEVKPNSSLERYDGIGPPFSCIFRIQNLGLFPIHGMMMKITIPIATRSGNRLLKLR DFLTDEANTSCNIWGNSTEYRPTPVEEDLRRAPQLNHSNSDVVSINCNIRLVPNQEIFHLLGNLWLRSLKALKYKSMKIMV NAALQRQFHSPFIFREEDPSRQIVFEISKQEDWQVP

靶点信息

研究背景	This gene encodes an alpha integrin. Integrins are heterodimeric integral membrane proteins composed of an alpha chain and a beta chain. This protein contains an I domain, is expressed in muscle tissue, dimerizes with beta 1 integrin in vitro, and appears to bind collagen in this form. Therefore, the protein may be involved in attaching muscle tissue to the extracellular matrix. Alternative transcriptional splice variants have been found for this gene, but their biological validity is not determined.
------	---

基因ID	22801
基因名	ITGA11
Swiss	Q9UKX5
别名	ITGA11;HsT18964

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn