

MAFB Rabbit pAb

货号: B20223

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	35kDa
实测分子量	50kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse liver,Mouse kidney,Rat liver
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 184-323 of human MAF B (NP_005452.2).
序列	PTSHPGPGPHATASATAAGGNGSVEDRFSDQLVSMVRELNRLRGFTKDEVIRLKQKRRTLKNRGYAQSCRYKRVQ QKHLENKTLIQQVEQLKQEVSRRLARERDAYKVKCEKLANSGFREAGSTSDSPSSPEFFL

靶点信息

研究背景	The protein encoded by this gene is a basic leucine zipper (bZIP) transcription factor that plays an important role in the regulation of lineage-specific hematopoiesis. The encoded nuclear protein represses ETS1-mediated transcription of erythroid-specific genes in myeloid cells. This gene contains no introns.
基因ID	9935

基因名	MAFB
Swiss	Q9Y5Q3
别名	MAFB;DURS3;KRML;MCTO

产品验证



Western blot analysis of MAFB expressed in Mouse liver,Mouse kidney,Rat liver using MAFB Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn