

IL17RD Rabbit pAb

货号: B20241

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	66kDa/79kDa/82kDa
实测分子量	110kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,HT-29,U-87MG,Mouse brain,Mouse kidney,Mouse lung,Rat brain
细胞定位	Cell membrane,Cytoplasm,Golgi apparatus membrane,Single-pass type I membrane protein,Single-pass t ype I membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 150-300 of human IL17 RD (NP_060033.3).
序列	FLNMKFETDYFVKVVPFPSIKNESNYHPFFFRTRACDLLLQPDNLACKPFWKPRNLNISQHGSDMQVSFDHAPHNFGFRF FYLHYKLKHEGPFKRKTCKQEQTETTSCLLQNVSPGDYIIELVDDTNTTRKVMHYALKPVHSPWAGPIRA

靶点信息

研究背景	This gene encodes a membrane protein belonging to the interleukin-17 receptor (IL-17R) protein family. T he encoded protein is a component of the interleukin-17 receptor signaling complex, and the interaction between this protein and IL-17R does not require the interleukin. The gene product also affects fibroblast growth factor signaling, inhibiting or stimulating growth through MAPK/ERK signaling. Alternate splicing g enerates multiple transcript variants encoding distinct isoforms.
------	--

基因ID	54756
基因名	IL17RD
Swiss	Q8NFM7
别名	IL17RD;HH18;IL-17RD;IL17RLM;SEF

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn