

MRPL32 Rabbit pAb

货号: B20250

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	21kDa
实测分子量	14kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,Raji,MCF7,HepG2,A-549,Mouse heart,Mouse kidney,Mouse liver
细胞定位	Mitochondrion
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 19-188 of human MRPL32 (NP_114109.1).
序列	VLRNYWERLLRKLPQSRPGFSPWPWGPALAVQGPAMFTEPANDTSGSKENSSLLDSIFWMAAPKNRRTIEVNRRCRRRNPPQ KLIKVKNNIDVCPEGHLKQKHVLCAYCEKVCKETAIEIRRIQIGKQEGGPFKAPTIEVTVVLYTGETPSEQDQGKRIIERDRKR PSWFTQN

靶点信息

研究背景	Mammalian mitochondrial ribosomal proteins are encoded by nuclear genes and help in protein synthesis within the mitochondrion. Mitochondrial ribosomes (mitoribosomes) consist of a small 28S subunit and a large 39S subunit. They have an estimated 75% protein to rRNA composition compared to prokaryotic ribosomes, where this ratio is reversed. Another difference between mammalian mitoribosomes and prokaryotic ribosomes is that the latter contain a 5S rRNA. Among different species, the proteins comprising the mitoribosome differ greatly in sequence, and sometimes in biochemical properties, which prevents easy recognition by sequence homology. This gene encodes a 39S subunit protein that belongs to the L32 ribosomal protein family. A pseudogene corresponding to this gene is found on chromosome Xp.
基因ID	64983
基因名	MRPL32
Swiss	Q9BYC8
别名	MRPL32;HSPC283;L32mt;MRP-L32;bMRP-59b

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn