

CD300A Rabbit pAb

货号: B20255

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:1000 - 1:5000
理论分子量	11kDa/16kDa/20kDa/33kDa
实测分子量	40-70KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	THP-1
细胞定位	Cell membrane,Single-pass type I membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 18-180 of human CD300 A (NP_009192.2).
序列	LSKCRTVAGPVGGSLSVQCPYEKEHRTLNKYWCRPPQIFLCDKIVETKGSAGKRNGRVSIRDSPANLSFTVTLENLTEEDA GTYWCGVDTPWLRDFHDPVVEVEVSVFPASTSMTPASITAAKTSTITTAFFPVSSSTTLFAVGATHSASIQEETEEVVNSQLP

靶点信息

研究背景	This gene encodes a member of the CD300 glycoprotein family of cell surface proteins found on leukocytes involved in immune response signaling pathways. This gene is located on chromosome 17 in a cluster with all but one of the other family members. Multiple transcript variants encoding different isoforms have been found for this gene.
基因ID	11314

基因名	CD300A
Swiss	Q9UGN4
别名	CD300A;CLM-8;CMRF-35-H9;CMRF-35H;CMRF35-H;CMRF35-H9;CMRF35H;CMRF35H9;IGSF12;IRC1;IRC1/IRC2;IRC2;IRp60

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn