

MYLK4 Rabbit pAb

货号: B20256

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	44kDa
实测分子量	40kDa/44kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,THP-1,U-937,HepG2,Mouse liver,Rat pancreas
细胞定位	
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 119-388 of human MYLK 4 (NP_001012418.2)
序列	QVHKCEETATGLKLAALKTRGMKDKEEVKNEISVMNQLDHANLIQLYDAFESKNDIVLVMMEYVDGGELFDRIIDESYNLTE LDTILFMKQICEGIRHMHQMYILHLDLKPENILCVNRDAKQIKIIDFGLARRYKPREKLKVNFGTPEFLAPEVVNYDFVSFPTD MWSVGVIAYMLLSGLSPFLGDNDAETLNNILACRWDLEDEEFQDISEEAKFISKLLIKEKSWRISASEALKHPWLSDHKLH SRLNAQKKKNRGSDAQDFVTK

靶点信息

研究背景	Predicted to enable myosin light chain kinase activity. Predicted to be involved in protein phosphorylation . Predicted to be integral component of membrane.
基因ID	340156

基因名	MYLK4
Swiss	Q86YV6
别名	MYLK4;SgK085

产品验证



Western blot analysis of MYLK4 expressed in Mouse liver,Rat pancreas,HeLa using MYLK4 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn