

Ezrin Rabbit pAb

货号: B20280

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC IP
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:100 IP: 1:50 - 1:100
理论分子量	69kDa
实测分子量	90kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	M21,SW480,U-251MG,293T,Mouse intestine,Mouse brain,Rat brain
细胞定位	Apical cell membrane,Cell projection,Cytoplasm,Cytoplasmic side,Peripheral membrane protein,cell cortex,cytoskeleton,microvillus membrane,ruffle membrane
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 350-470 of human Ezrin in (NP_003370.2).
序列	RLQDYEEKTKKAERELSEQIQRALQLEEERKRAQEEAERLEADRMAALRAKEELERQAVDQIKSSEQLAAELA EYTAKIALLE EARRRKEDEVEEWQHRAKEAQDDL VKTKEELHLVMTAP

靶点信息

研究背景	The cytoplasmic peripheral membrane protein encoded by this gene functions as a protein-tyrosine kinase substrate in microvilli. As a member of the ERM protein family, this protein serves as an intermediate between the plasma membrane and the actin cytoskeleton. This protein plays a key role in cell surface structure adhesion, migration and organization, and it has been implicated in various human cancers. A pseudogene located on chromosome 3 has been identified for this gene. Alternatively spliced variants have also been described for this gene.
基因ID	7430
基因名	EZR
Swiss	P15311
别名	CVIL;CVL;HEL-S-105;VIL2;EZR;Ezrin;ezrin

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn