

TNFRSF6B Rabbit pAb

货号: B20286

产品信息

反应	Human,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	32kDa
实测分子量	33KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,A-549,DU145,SGC-7901,SW480
细胞定位	Secreted
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 200-300 of human TNFRSF6B (NP_003814.1).
序列	PLSTRVPGAEECERAVIDFVAFQDISIKRLQRLLQALEAPEGWGPTPRAGRAALQLKLRRRLTELLGAQDGALLVRLLQALRVARMPGLERSVRERFLPVH

靶点信息

研究背景	This gene belongs to the tumor necrosis factor receptor superfamily. The encoded protein is postulated to play a regulatory role in suppressing FasL- and LIGHT-mediated cell death. It acts as a decoy receptor that competes with death receptors for ligand binding. Over-expression of this gene has been noted in gastrointestinal tract tumors. Read-through transcription into this gene from the neighboring upstream gene, which encodes regulator of telomere elongation helicase 1 (RTEL1), generates a non-coding transcript.
基因ID	8771

基因名	TNFRSF6B
Swiss	O95407
别名	TNFRSF6B;DCR3;DJ583P15.1.1;M68;M68E;TR6

产品验证



Western blot analysis of TNFRSF6B expressed in HepG2,A-549,DU145 using TNFRSF6B Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn