

DOK4 Rabbit pAb

货号: B20323

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	37kDa
实测分子量	37KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse heart,Mouse kidney,Rat heart,Rat kidney
细胞定位	cytoplasm,cytosol
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-326 of human DOK4 (NP_060580.2).
序列	MATNFSDIVKQGYVKMKSRKLGIIYRRCWLVFRKSSSKGPQRLEKYPDEKSVCLRGCPKVTEISNVKCVTRLPKETKRQAVA IIFTDDSARTFTCDSELEAEWYKTLSEVCLGSRLNDISLGEPDLLAPGVQCEQTDRFNVFLLPCPNLDVYGECKLQITHENI YLWDIHNPVKLVSWPLCSLRRYGRDATRFTFEAGRMCDAGEGLYTFQTQEGEQIYQRVHSATLAI AEQH KRVLLEMEKN VRLLNKGTEHYSYPCTPTTMLPR SAYWHHITGSQNI AEASSYAGEGYGAAQASSETDLLNRFILLKPKPSQGD SSEAKTPS Q

靶点信息

研究背景	DOK proteins are enzymatically inert adaptor or scaffolding proteins. They provide a docking platform for the assembly of multimolecular signaling complexes. DOK4 functions in RET-mediated neurite outgrowth and plays a positive role in activation of the MAP kinase pathway (By similarity). Putative link with downstream effectors of RET in neuronal differentiation. May be involved in the regulation of the immune response induced by T-cells.
基因ID	55715
基因名	DOK4
Swiss	Q8TEW6
别名	DOK4;IRS-5;IRS5

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn