

# Phospho-EIF4G1-S1108 Rabbit pAb

货号: B20386

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	154kDa/158kDa/166kDa/171kDa/175kDa/176kDa
实测分子量	250kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	NIH/3T3,NIH/3T3+Serum
细胞定位	cytoplasm,cytoplasmic stress granule,cytosol,eukaryotic translation initiation factor 4F complex,nucleus, postsynapse
纯化	Affinity purification

抗原信息


抗原信息	A synthetic phosphorylated peptide around S1108 of human EIF4G1 (NP_001278086.1).
序列	SLSRE

靶点信息

研究背景	The protein encoded by this gene is a component of the multi-subunit protein complex EIF4F. This complex facilitates the recruitment of mRNA to the ribosome, which is a rate-limiting step during the initiation phase of protein synthesis. The recognition of the mRNA cap and the ATP-dependent unwinding of 5'-terminal secondary structure is catalyzed by factors in this complex. The subunit encoded by this gene is a large scaffolding protein that contains binding sites for other members of the EIF4F complex. A domain at its N-terminus can also interact with the poly(A)-binding protein, which may mediate the circularization of mRNA during translation. Alternative splicing results in multiple transcript variants, some of which are derived from alternative promoter usage.
------	---

基因ID	1981
基因名	EIF4G1
Swiss	Q04637
别名	EIF4G1;EIF-4G1;EIF4F;EIF4G;EIF4GI;P220;PARK18

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)