

DPPA5 Rabbit pAb

货号: B20435

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	13kDa
实测分子量	13kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse heart,Rat testis
细胞定位	cytoplasm
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-116 of human DPPA5 (NP_001020461.1).
序列	MGTLPARRHIPPWVKVPEDLKDPEVFQVQTRLLKAIFGPDGSRIPYIEQVSKAMLELKALESSDLTEVVVYGSYLYKLRTKW MLQSMAEWHRRQRQERGMLKLAEAMNALELGPWMK

靶点信息

研究背景	This gene encodes a protein that may function in the control of cell pluripotency and early embryogenesis . Expression of this gene is a specific marker for pluripotent stem cells. Pseudogenes of this gene are located on the short arm of chromosome 10 and the long arm of chromosomes 14 and 19.
基因ID	340168

基因名	DPPA5
Swiss	A6NC42
别名	DPPA5;ESG1

产品验证



Western blot analysis of DPPA5 expressed in Mouse heart,Rat testis using DPPA5 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn