

MYCL Rabbit pAb

货号: B20448

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	24kDa/40kDa/43kDa
实测分子量	40KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,BxPC-3
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 40-180 of human MYCL (NP_005367.2).
序列	FYDYDCGEDFYRSTAPSEDIWKKFELVPSPTSPPWGLGPGAGDPAPGIGPPEPWGGCTGDEAESRGHSGWGRNYA SIIRRDCMWSGFSARERLERAVSDRLAPGAPRGNPPKASAAPDCTPSLEAGNPAPAAPCPLGE

靶点信息

研究背景	Predicted to enable DNA-binding transcription factor activity, RNA polymerase II-specific and RNA polyme rase II cis-regulatory region sequence-specific DNA binding activity. Predicted to be involved in regulation of transcription by RNA polymerase II. Predicted to act upstream of or within regulation of inner ear audit ory receptor cell differentiation. Located in chromosome and nucleoplasm.
基因ID	4610

基因名	MYCL
Swiss	P12524
别名	L-Myc;LMYC;MYCL1;bHLHe38;MYCL

产品验证



Western blot analysis of MYCL expressed in A-549,BxPC-3 using MYCL Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn