

SON Rabbit pAb

货号: B20454

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	264kDa
实测分子量	250kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	LO2,HeLa
细胞定位	nuclear speck
纯化	Affinity purification

抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 850-950 of human SON (NP_001278 340.1).
序列	TMDSQMLATSSMDSQMLASGTMDSQMLASGTMDAQMLASGTMDAQMLASSTQDSAMLGSKSPDPYRLAQDPYRLAQDPYRLGHD PYRLGHDAYRLGQDPYR

靶点信息

研究背景	This gene encodes a protein that contains multiple simple repeats. The encoded protein binds RNA and promotes pre-mRNA splicing, particularly of transcripts with poor splice sites. The protein also recognizes a specific DNA sequence found in the human hepatitis B virus (HBV) and represses HBV core promoter activity. There is a pseudogene for this gene on chromosome 1. Alternative splicing results in multiple transcript variants.
基因ID	6651

基因名	SON
Swiss	P18583
别名	SON;BASS1;C21orf50;DBP-5;NREBP;SON3;TOKIMS

产品验证



Western blot analysis of SON expressed in LO2,HeLa using SON Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn