

DPF3 Rabbit pAb

货号: B20483

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	43kDa
实测分子量	43KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse heart
细胞定位	nucleoplasm
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 90-200 of human DPF3 (NP_036206.3).
序列	DPKLRLLLEIKPEVELPLKKGFTSESTTLEALLRGEGVEKKVDAREEESIQEIQRVLENDENVEEGNEEEDLEEDIPKRKNRTRGRARGSAGGRRRRHDAASQEDHDKPYVC

靶点信息

研究背景	This gene encodes a member of the D4 protein family. The encoded protein is a transcription regulator th at binds acetylated histones and is a component of the BAF chromatin remodeling complex. Alternate spli cing results in multiple transcript variants encoding different isoforms.
基因ID	8110

基因名	DPF3
Swiss	Q92784
别名	DPF3;BAF45C;CERD4

产品验证



Western blot analysis of DPF3 expressed in Mouse heart using DPF3 Rabbit pAb at 1:1000. Secondary anti body: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% no n-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn