

GNA15 Rabbit pAb

货号: B20546

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	43kDa
实测分子量	43kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	LO2,HeLa,mouse spleen
细胞定位	plasma membrane,synapse
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 112-374 of human GNA 15 (NP_002059.3).
序列	HASLVMSQDPYKVTTFEKRYAAAMQWLWRDAGIRACYERRREFHLLDSAVYYLSHLERITEEGYVPTAQDVLRSRMPPTTGI NEYCFSVQKTNLRIVDVGGQKSERKKWIHCFENVIALIYLASLSEYDQCLEENNQENRMKESLALFGTILELPWFKSTSVILF LNKTDILEEKIPTSHLATYFPSFQGPQDAEAAKRFILDMYTRMYTGCVDGPEGSKKGARSRRLLFSHYTCATDTQNIRKVFK DVRDSVLARYLDEINLL

靶点信息

研究背景	Guanine nucleotide-binding proteins (G proteins are involved as modulators or transducers in various tran smembrane signaling systems.
基因ID	2769

基因名	GNA15
Swiss	P30679
别名	GNA15;GNA16

产品验证



Western blot analysis of GNA15 expressed in LO2,HeLa,mouse spleen using GNA15 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn