

IFNA21 Rabbit pAb

货号: B20555

产品信息

反应	Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	
实测分子量	22kD
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	extracellular region,extracellular space
纯化	Affinity purification

抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 100 to the C-terminus of human IFN A21 (NP_002166.2).
序列	WEQSLLEKFSTELNQQLNDLEACVIQEVGVEETPLMNVDSILAVKKYFQRITLYLTEKKYSPCAWEVVRAEIMRSFSLSKIFQERLRRKE

靶点信息

研究背景	This gene is a member of the alpha interferon gene cluster on the short arm of chromosome 9. Interferon s are cytokines produced in response to viral infection that mediate the immune response and interfere w ith viral replication. The encoded protein is a type I interferon and may play a specific role in the antiviral response to rubella virus.
基因ID	3452

基因名	IFNA21
Swiss	P01568
别名	IFNA21;IFN-alpha1;LeIF F;leIF-F

产品验证



Western blot analysis of IFNA21 expressed in mouse liver using IFNA21 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn