

THBS1 Rabbit pAb

货号: B20607

产品信息

反应	Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	IHC
推荐浓度	IHC: 1:50 - 1:200
理论分子量	120kDa/129kDa
实测分子量	120kDa/129kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Endoplasmic reticulum,Sarcoplasmic reticulum
纯化	Affinity purification

抗原信息

抗原信息	A synthetic peptide corresponding to a sequence within amino acids 150-250 of human THBS1 (NP_003237.2).
序列	TGQWKSITLFVQEDRAQLYIDCEKMENAELDVPIQSVFTRDLASIALRLIAKGGVNDNFQGV LQNVRFVFGTTPEDILRNKGCSSSTSVLLTLDNNVVNGS

靶点信息

研究背景	The protein encoded by this gene is a subunit of a disulfide-linked homotrimeric protein. This protein is an adhesive glycoprotein that mediates cell-to-cell and cell-to-matrix interactions. This protein can bind to fibrinogen, fibronectin, laminin, type V collagen and integrins alpha-V/beta-1. This protein has been shown to play roles in platelet aggregation, angiogenesis, and tumorigenesis.
基因ID	7057

基因名	THBS1
Swiss	P07996
别名	THBS1;THBS;THBS-1;TSP;TSP-1;TSP1

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn