

SCGB1A1 Rabbit pAb

货号: B20612

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	9kDa
实测分子量	11kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse lung
细胞定位	Secreted
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 22-91 of human SCGB1 A1 (NP_003348.1).
序列	EICPSFQRVIETLLMDTPSSYEAAAMELFSPDQDMREAGQLKKLVDTLPQKPRESIIKLMEKIAQSSLCN

靶点信息

研究背景	This gene encodes a member of the secretoglobin family of small secreted proteins. The encoded protein has been implicated in numerous functions including anti-inflammation, inhibition of phospholipase A2 and the sequestering of hydrophobic ligands. Defects in this gene are associated with a susceptibility to asthma.
基因ID	7356

基因名	SCGB1A1
Swiss	P11684
别名	SCGB1A1;CC10;CC16;CCPBP;CCSP;UGB;UP-1;UP1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn