

ZNF711 Rabbit pAb

货号: B20618

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	
实测分子量	110kD
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	nucleus
纯化	Affinity purification

抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 200-300 of human ZNF711 (NP_001317503.1).
序列	SEDYLMISLDDVGEKLEHMGNTPLKIGSDGSQEDAKEDGFGSEVIKVYIFKAEAEDDVEIGGTEIVTESEYTSGHSVAGVL DQSRMQREKMVYMAVKDSSQ

靶点信息

研究背景	This gene encodes a zinc finger protein of unknown function. It bears similarity to a zinc finger protein which acts as a transcriptional activator. This gene lies in a region of the X chromosome which has been associated with mental retardation.
基因ID	7552

基因名	ZNF711
Swiss	Q9Y462
别名	ZNF711;CMPX1;MRX97;ZNF4;ZNF5;ZNF6;Zfp711;dj75N13.1

产品验证



Western blot analysis of ZNF711 expressed in human breast cancer using ZNF711 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn