

WNT7B Rabbit pAb

货号: B20623

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	39kDa
实测分子量	39kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,HepG2
细胞定位	
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 120-349 of human WNT 7B (NP_478679.1).
序列	TAACSQGNLSNCGCDREKQGYYNQAEGWKWGGCSADVRYGIDFSRRFVDAREIKKNARRLMNLHNNEAGRKVLEDRLM QLECKCHGVSGSCTTKTCWTTLPKFREVGHLLKEKYNAAVQVEVVRASRLRQPTFLRIKQLRSYQKPMETDLVYIEKSPNY CEEDAATGSVGTTQGRLCNRTSPGADGCDTMCCGRGYNTHQYTKVWQCNCCKFWCCFVKCNTCSERVEFTCK

靶点信息

研究背景	This gene is a member of the WNT gene family, which consists of structurally related genes that encode s ecreted signaling proteins. These proteins have been implicated in oncogenesis and in several developme ntal processes, including regulation of cell fate and patterning during embryogenesis. Among members of the human WNT family, this gene product is most similar to WNT7A protein.
基因ID	7477

基因名	WNT7B
Swiss	P56706
别名	WNT7B

产品验证



Western blot analysis of WNT7B expressed in A-549,HepG2 using WNT7B Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn