

# MAML1 Rabbit pAb

货号: B20644

产品信息

|       |                                                                                                      |
|-------|------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse,Rat                                                                                      |
| 宿主    | Rabbit                                                                                               |
| 克隆性   | Polyclonal                                                                                           |
| 预测反应  |                                                                                                      |
| 应用    | WB                                                                                                   |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000                                                                            |
| 理论分子量 | 108kDa                                                                                               |
| 实测分子量 | 120kDa                                                                                               |
| 形式    | Liquid                                                                                               |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                         |
| 阳性对照  | HeLa,Jurkat,mouse brain,rat brain                                                                    |
| 细胞定位  | nuclear speck,nucleoplasm,nucleus                                                                    |
| 纯化    | Affinity purification                                                                                |

抗原信息


|      |                                                                                                                                                                                                                                                                                                                       |
|------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human MAML1 (NP_055572.1).                                                                                                                                                                                                     |
| 序列   | MVLPTCPMAEFALPRHSVMERLRRRIELCRRHHSTCEARYEAVSPERLELERQHTFALHQRCIQAKAKRAGKHRQPPAATA<br>PAPAAPAPRLDAADGPEHGRPATHLHDTVKNLDSATSPQNGDQQNGYGDLFPGHKKTRREAPLGVAISSNGLPPASPL<br>GQSDKPSGADALQSSGKHSGLDSLNNKKRLADSSLHLNGGSGNPSEFPLSLNKEKQEPVEDLPCMITGTGVSISQSNLM<br>PDLNLNEQEWKELIEELNRSVPDEDMKDLFNEDFEEKDPESGSGSATQTPLAQDINIKT |

靶点信息

|      |                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
|------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | This protein is the human homolog of mastermind, a Drosophila protein that plays a role in the Notch signaling pathway involved in cell-fate determination. There is in vitro evidence that the human homolog forms a complex with the intracellular portion of human Notch receptors and can increase expression of a Notch-induced gene. This evidence supports its proposed function as a transcriptional co-activator in the Notch signaling pathway. |
|------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|       |                  |
|-------|------------------|
| 基因ID  | 9794             |
| 基因名   | MAML1            |
| Swiss | Q92585           |
| 别名    | MAML1;Mam-1;Mam1 |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)