

PRPF6 Rabbit pAb

货号: B20672

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	102kDa/106kDa
实测分子量	102kDa/106kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 592-941 of human PRPF6 (NP_036601.2).
序列	SLEALLQRAVAHCPKAEVLWLMGAKSKWLAGDVPAAARSILALAFQANPNSEIWLAAVKLESENDEYERARRLLAKARSSAPTARVFMKSVKLEWVQDNIRAAQDLCEEALRHYEDFPKLWMMKGQIEEQKEMMEKAREAYNQGLKKCPHSTPLWLLLSRLEEKIGQLTRARAILEKSRLKNPKNPGLWLESVRLEYRAGLKNIANTLMAKALQECPSGILWSEAIFLEARPQRRTKSVDALKKCEHDPHVLLAVAKLFWSQLRKITKAREWFHRTVKIDSDLGDAWAFFYKFELQHGTEEQQEEVRKRCESAEPRHGELWC AVSKDIANWQKKIGDILRLVAGRIKNTF

靶点信息

研究背景	The protein encoded by this gene appears to be involved in pre-mRNA splicing, possibly acting as a bridging factor between U5 and U4/U6 snRNPs in formation of the spliceosome. The encoded protein also can bind androgen receptor, providing a link between transcriptional activation and splicing.
------	--

基因ID	24148
基因名	PRPF6
Swiss	O94906
别名	PRPF6;ANT-1;ANT1;C20orf14;Prp6;RP60;SNRNP102;TOM;U5-102K;hPrp6

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn