

CNNM1 Rabbit pAb

货号: B20675

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:100
理论分子量	
实测分子量	104kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain,Rat brain
细胞定位	plasma membrane
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 692-951 of human CNN M1 (NP_065081.2).
序列	NGAFTYYGVPAIMTTACSDNDVRKVGSLAGSSVFLNRSPSRCSGLNRSESPNRERSDFGGSNTQLYSSSNLYMPDYSV HILSDVQFVKITRQQYQNALTACHMDSSPQSPDMEAFTDGDSTKAPTTRGTPQTPKDDPAITLLNNRNSLPCSRSDGLRS PSEVVYLRMEELAFTQEEMTDFEEHSTQQLTSPAAPVTRAASDSECCNINLDTETSPCSSDFEENVGKKLLRTLSGQKRK RSPGERTSEDNSNLTPLIT

靶点信息

研究背景	This gene encodes a member of the ancient conserved domain protein family. The encoded protein may b ind copper. Alternative splicing results in multiple transcript variants.
基因ID	26507

基因名	CNNM1
Swiss	Q9NRU3
别名	CNNM1;ACDP1;CLP-1

产品验证



Western blot analysis of CNNM1 expressed in Mouse brain,Rat brain using CNNM1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn