

CAMKV Rabbit pAb

货号: B20724

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	51kDa/54kDa
实测分子量	54kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse eye,Mouse spleen,Rat eye
细胞定位	glutamatergic synapse,plasma membrane
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-285 of human CAMKV (NP_076951.2).
序列	MPFGCVTLGDKKNYNPSEVTDRLDGLQVIKTEEFCEIFRAKDKTTGKLHTCKKFQKRDGRKVRKAAKNEIGILKMKHPN ILQLVDVFTVTRKEYFIFLELATGREVFDWILDQGYYSERDTSNVVRQVLEAVAYLHSLKIVHRNLKLENLVYYNRLKNSKIVIS DFHLAKLENGLIKEPCGTPPEYLAPEVVGRQRYGRPVDCWAIGVIMYILLSGNPPFYEEVEEDDYENHDKNLFRKILAGDYEF DSPYWDDISQAADLVTRLMEVEQDQRITAEAAISHEW

靶点信息

研究背景	Does not appear to have detectable kinase activity.
基因ID	79012

基因名	CAMKV
Swiss	Q8NCB2
别名	CAMKV;1G5;VACAMKL

产品验证



Western blot analysis of CAMKV expressed in Mouse eye,Mouse spleen,Rat eye using CAMKV Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn