

PPP1R9B Rabbit pAb

货号: B20734

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	89kDa
实测分子量	127kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,K-562,Mouse brain
细胞定位	Cell junction,Cell membrane,Cell projection,Cytoplasm,Nucleus,adherens junction,cytoskeleton,dendritic spine,filopodium,lamellipodium,ruffle membrane,synapse
纯化	Affinity purification

抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 250-350 of human PPP1R9B (NP_115984.3).
序列	VFQPPPPPPAPSGDAPAEKERCPAGQQPPQHRVAPARPPPKPREVRKIKPVEVEESGESEAESAPGEVIQAEVTVHAALEN GSTVATAASPAPEEPAQA

靶点信息

研究背景	This gene encodes a scaffold protein that functions as a regulatory subunit of protein phosphatase 1a. Expression of this gene is particularly high in dendritic spines, suggesting that the encoded protein may play a role in receiving signals from the central nervous system. The encoded protein has putative tumor suppressor function and decreased expression has been observed in tumors.
------	---

基因ID	84687
基因名	PPP1R9B
Swiss	Q96SB3
别名	PPP1R9B;PPP1R6;PPP1R9;SPINO;Spn

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn