

# TRIP13 Rabbit pAb

货号: B20784

## 产品信息

|       |                                                                                                      |
|-------|------------------------------------------------------------------------------------------------------|
| 反应    | Human                                                                                                |
| 宿主    | Rabbit                                                                                               |
| 克隆性   | Polyclonal                                                                                           |
| 预测反应  |                                                                                                      |
| 应用    | WB                                                                                                   |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000                                                                            |
| 理论分子量 | 22kDa/48kDa                                                                                          |
| 实测分子量 | 48kDa                                                                                                |
| 形式    | Liquid                                                                                               |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                         |
| 阳性对照  | HeLa,PC-3                                                                                            |
| 细胞定位  | male germ cell nucleus,nucleus                                                                       |
| 纯化    | Affinity purification                                                                                |

## 抗原信息


|      |                                                                                                                            |
|------|----------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human TRIP13 (NP_004228.1).         |
| 序列   | MDEAVGDLKQALPCVAESPTVHVEVHQRGSSSTAKKEDINLSVRKLLNRHNIVFGDYTWTEFDEPFLTRNVQSVSIIDTELKVKDSQPIDLSACTVALHIFQLNEDGPSSSENLEEETENII |

## 靶点信息

|      |                                                                                                                                                                                                                                                                                                       |
|------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | This gene encodes a protein that interacts with thyroid hormone receptors, also known as hormone-dependent transcription factors. The gene product interacts specifically with the ligand binding domain. This gene is one of several that may play a role in early-stage non-small cell lung cancer. |
| 基因ID | 9319                                                                                                                                                                                                                                                                                                  |

|       |               |
|-------|---------------|
| 基因名   | TRIP13        |
| Swiss | Q15645        |
| 别名    | TRIP13;16E1BP |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)