

FILIP1L Rabbit pAb

货号: B20824

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	130kDa
实测分子量	110kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HT-29,mouse ovary,A375
细胞定位	cytoplasm,nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human FILIP1L (NP_878913.2).
序列	MRSRGSDTEGSAQKKFPRHTKGHSFQGPKNMKHRQQDKDSPESDVILPCPKAEKPHSGNGHQAEDLSRDDLLFLLSILEGELQARDEVIGILKAEKMDLALLEAQYGFVTPKKVLEALQRDAFAQSTPWQEDIYEKPMNELDKVVEKHKESYRRILGQLLVAEKSRRQTILELEEEK RKHKEYMEKSDEFICLLEQEC

靶点信息

研究背景	Acts as a regulator of the antiangiogenic activity on endothelial cells. When overexpressed in endothelial cells, leads to inhibition of cell proliferation and migration and an increase in apoptosis. Inhibits melanom a growth When expressed in tumor-associated vasculature.
基因ID	11259

基因名	FILIP1L
Swiss	Q4L180
别名	DOC1;DOC-1;GIP90;GIP130;FILIP1L

产品验证



Western blot analysis of FILIP1L expressed in HT-29,mouse ovary,A375 using FILIP1L Rabbit pAb at 1:10 00. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Bl ocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn