

SV2A Rabbit pAb

货号: B20826

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	83kDa
实测分子量	105kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse kidney,Mouse heart,Mouse brain,Rat brain
细胞定位	dendrite,endoplasmic reticulum,GABA-ergic synapse,glutamatergic synapse,neuromuscular junction,neur on projection,plasma membrane,presynaptic active zone,synaptic vesicle,synaptic vesicle membrane
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 356-447 of human SV2A (NP_055664.3).
序列	ESPRFFLENGKHDEAWMLKQVHDTNMRAKGHPERVFSVTHIKTIHQEDELIEIQSDTGTWYQRWGVRLSLGGQVWG NFLSCFGPEYRRIT

靶点信息

研究背景	The protein encoded by this gene is one of three related synaptic vesicle proteins. The encoded protein may interact with synaptotagmin to enhance low frequency neurotransmission in quiescent neurons.
基因ID	9900

基因名	SV2A
Swiss	Q7L0J3
别名	SV2;SV2A

产品验证



Western blot analysis of SV2A expressed in Mouse kidney,Rat brain,Mouse heart using SV2A Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn