

SIGLEC1 Rabbit pAb

货号: B20838

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	200kDa
实测分子量	200kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2
细胞定位	early endosome,extracellular region,late endosome,plasma membrane
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 50-300 of human SIGLE C1 (NP_075556.1).
序列	EVPDGTAIWYYDYSGQRQVVSADPKLVEARFRGRTEFMGNPEHRVCNLLKDLQPEDSGSYNFRFEISEVNRWSDVK GTLVTVTEEPVPTIASPVELLEGTEVDFNCSTPYVCLQEQVRLQWQGQDPARSVTFNSQKFEPTGVGHLETLHMAMSW QDHGRILRCQLSVANHRAQSEIHLQVKYAPKGVKILLSPSGRNILPGELVLTLCQVNSSYPAVSSIKWLKDGVRQLQTKTGVL HLPQAAWSDA

靶点信息

研究背景	This gene encodes a member of the immunoglobulin superfamily. The encoded protein is a lectin-like adhesion molecule that binds glycoconjugate ligands on cell surfaces in a sialic acid-dependent manner. It is a type I transmembrane protein expressed only by a subpopulation of macrophages and is involved in mediating cell-cell interactions. Alternative splicing produces a transcript variant encoding an isoform that is soluble rather than membrane-bound; however, the full-length nature of this variant has not been determined.
------	---

基因ID	6614
基因名	SIGLEC1
Swiss	Q9BZZ2
别名	SN;CD169;SIGLEC-1;SIGLEC1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn