

HEYL Rabbit pAb

货号: B20840

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	328kDa
实测分子量	37KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U266,MDA-MB-231,Jurkat,Molt-4(Negative control)
细胞定位	cytoplasm,nucleoplasm,nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 129-328 of human HEYL (NP_055386.1).
序列	IRYLGVLGEPSSRADPVRI RL LSHLNSYAAEMEPSPTPTG PLAFPWPWSFFHSCPLPALSNQLAILGRVPSPVLPGVSSP AYIPALRTAPLRRATGIILPARRNVLP SRGASSTRRARPLER PATPVVPAPSSRAARSSHIAPLLQSSSPTPPGPTGSAAYVA VPTPNSSSPGPAGRPAGAMLYHSWVSEITEIGAF

靶点信息

研究背景	This gene encodes a member of the hairy and enhancer of split-related (HESR) family of basic helix-loop-helix (bHLH)-type transcription factors. The sequence of the encoded protein contains a conserved bHLH and orange domain, but its YRPW motif has diverged from other HESR family members. It is thought to be an effector of Notch signaling and a regulator of cell fate decisions. Alternatively spliced transcript variants have been found, but their biological validity has not been determined.
------	--

基因ID	26508
基因名	HEYL
Swiss	Q9NQ87
别名	HEY3;HRT3;HESR3;bHLHb33;HEYL

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn