

DNAH5 Rabbit pAb

货号: B20867

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	
实测分子量	529KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse testis
细胞定位	9+0 motile cilium,9+2 motile cilium,axoneme,cytoplasm,extracellular region
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 400-650 of human DNA H5 (NP_001360.1).
序列	VTNQIISACKAYITNNGTASIWNQPQDVVEEKILSAIKLKQEYQLCFHKTKQKLKQNPNAKQDFSEMYIFGKFETFHRRLLAKIIDIFTTLKTYSVLQDSTIEGLEDMATKYQGIVATIKKKEYNFLDQRKMDFDQDYEEFCKQTNDLHNELRKFMVDVTFAKIQNTNQALRMLKKFERLNIPNLGIDDKYQLILENYGADIDMISKLYTKQKYPPLARNQPPIAGKILWARQLFHRIQQPMLFQQHPAV

靶点信息

研究背景	This gene encodes a dynein protein, which is part of a microtubule-associated motor protein complex consisting of heavy, light, and intermediate chains. This protein is an axonemal heavy chain dynein. It functions as a force-generating protein with ATPase activity, whereby the release of ADP is thought to produce the force-producing power stroke. Mutations in this gene cause primary ciliary dyskinesia type 3, as well as Kartagener syndrome, which are both diseases due to ciliary defects. [provided by RefSeq, Oct 2009]
------	---

基因ID	1767
基因名	DNAH5
Swiss	Q8TE73
别名	CILD3;DNAHC5;HL1;KTG NR;PCD;DNAH5

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn