

ITPKA Rabbit pAb

货号: B20888

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	51kDa
实测分子量	55kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain,Rat brain
细胞定位	cytoplasm,cytoskeleton,cytosol,dendritic spine,glutamatergic synapse,nucleus,postsynaptic actin cytoske leton
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 141-461 of human ITPK A (NP_002211.1).
序列	LLSDSESRSRGNVQLEAGEDVGQKNHWQKIRTMVNLPVISPFKKRYAWVQLAGHTGSFKAAGTSGLILKRCSEPERYCLA RLMADALRGCVPAFHGVVERDGESYLQLQDLLDGF DGPCVLDCKMGVRTYLEEELTKARERPKLRKDMYKKMLAVDPEA PTEEEHAQRAVTKPRYMQWREGISSSTTLGFRIEGIKKADGSCSTDFKTTRSREQVLRVFEEFVQGDDEEVLRRYLNRLQQI RDTLEVSEFFRRHEVIGSSLLFVHDHCHRAGVWLIDFGKTTPLPDGQILDHRRPWEEGNREDGYLLGLDNLIGILASLAER

靶点信息

研究背景	Regulates inositol phosphate metabolism by phosphorylation of second messenger inositol 1,4,5-trisphosphate to Ins(1,3,4,5)P4. The activity of the inositol 1,4,5-trisphosphate 3-kinase is responsible for regulating the levels of a large number of inositol polyphosphates that are important in cellular signaling. Both calcium/calmodulin and protein phosphorylation mechanisms control its activity. It is also a substrate for the cyclic AMP-dependent protein kinase, calcium/calmodulin-dependent protein kinase II, and protein kinase C in vitro.[provided by RefSeq, Apr 2011]
基因ID	3706
基因名	ITPKA
Swiss	P23677
别名	IP3-3KA;IP3KA;ITPKA

产品验证



Western blot analysis of ITPKA expressed in Mouse brain,Rat brain using ITPKA Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn