

SLC9A3R2 Rabbit pAb

货号: B20938

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	
实测分子量	50kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa
细胞定位	apical plasma membrane,extracellular exosome,focal adhesion,nucleus,plasma membrane
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 180-300 of human SLC9 A3R2 (NP_001123484).
序列	DPGSPAARSGRLRAQDRLLIEVNGQNVEGLRHAEEVVASIKAREDEARLLVDPETDEHFKRLRVTPTEEHVEGPLPSPVTNGT SPAQLNGGSACSSRSDLPGSDKDTEGSAWKQDPFQESGL

靶点信息

研究背景	This gene encodes a member of the NHERF family of PDZ scaffolding proteins. These proteins mediate many cellular processes by binding to and regulating the membrane expression and protein-protein interactions of membrane receptors and transport proteins. The encoded protein plays a role in intestinal sodium absorption by regulating the activity of the sodium/hydrogen exchanger 3, and may also regulate the cystic fibrosis transmembrane regulator (CFTR) ion channel. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Nov 2011]
------	--

基因ID	9351
基因名	SLC9A3R2
Swiss	Q15599
别名	E3KARP;NHE3RF2;NHERF-2;NHERF2;OCTS2;SIP-1;SIP1;TKA-1;SLC9A3R2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn