

RBM19 Rabbit pAb

货号: B20944

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	107kDa
实测分子量	107kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Raji
细胞定位	cytoplasm,nuclear speck,nucleolus,nucleoplasm
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 759-960 of human RBM 19 (NP_057280.2).
序列	CSISKKKNKAGVLLSMGFGFVEYRKPEQAQKALKQLQGHVVDGHKLEVRISERATKPAVTLARKKQVPRKQTTSKILVRNIP FQAHSREIRELFSTFGELKTVRLPKKMTGTGTHRGFGFVDFLTQDAKRAFNALCHSTHLYGRRLVLEWADSEVTLQALR RKTAAHFHEPPKKKRSVVLDEILEQLEGSDSDSEEQTLQL

靶点信息

研究背景	This gene encodes a nucleolar protein that contains six RNA-binding motifs. The encoded protein may be involved in regulating ribosome biogenesis. Multiple alternatively spliced variants, encoding the same protein, have been identified.[provided by RefSeq, Apr 2009]
基因ID	9904

基因名	RBM19
Swiss	Q9Y4C8
别名	RBM19

产品验证



Western blot analysis of RBM19 expressed in Raji using RBM19 Rabbit pAb at 1:1000. Secondary antibody : HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn