

SUPT7L Rabbit pAb

货号: B20945

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	30kDa/45kDa/46kDa
实测分子量	60kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Rat lung
细胞定位	nucleoplasm,nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 210-400 of human SUPT 7L (NP_001269658).
序列	TPFPDVMEQVFHEVGIGSVLSLQKFWQHRIKDYHSYMLQISKQLSEEYERIVNPEKATEDAKPVKIKEEPVSDITFPVSEELEADLASGDQSLPMGVLGAQSERFPSNLEVEASPQASSAEVNASPLWNLAHVKMEPQESEEENVS GHGVLGSDVFEEPMSGMSEAGIPQSPDDSDSSYGSHSTDSLMGSS

靶点信息

研究背景	SUPT7L is a protein subunit of the human STAGA complex (SPT3; (MIM 602947)/TAF9 (MIM 600822)/GCN5 (MIM 602301) acetyltransferase complex), which is a chromatin-modifying multiprotein complex (Martinez et al., 2001 [PubMed 11564863]).[supplied by OMIM, Apr 2009]
基因ID	9913

基因名	SUPT7L
Swiss	O94864
别名	SPT7L;STAF65;STAF65(gamma);STAF65G;SUPT7H;SUPT7L

产品验证



Western blot analysis of SUPT7L expressed in Rat lung using SUPT7L Rabbit pAb at 1:1000. Secondary anti body: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% no n-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn